

## SUCCESS STORIES

**Specific Genome Editing by TALEN and CRISPR-Cas9 technology** Transcription activator-like effector nucleases (TALENs) are artificial restriction enzymes generated by fusing a TAL effector DNA binding domain to a DNA cleavage domain. TAL Nuclease has revolutionized the array of genome editing based on its ability to create a site-specific double-stranded DNA break (DSB). The break is then subsequently repaired by cellular machinery, through either homology-dependent repairs or non-homologous end joining (NHEJ), that creates a deletion mutation or a knock-in site for introducing site-specific DNA sequences of interest. Clustered, regularly interspaced, short palindromic repeat (CRISPR) RNA-guided nucleases (RGN) can robustly induce genome editing. Repair of RGN-induced double-stranded breaks by non-homologous end-joining or homology repair introduces insertion or deletion mutations or specific sequence alterations. The *Streptococcus pyogenes* Cas9 nuclease (Cas9) cleaves the intervening spacer sequence directed by a single guide RNA (gRNA) with a 20 nucleotides (nt) target complementarity region at its 5' end.




Table 1 describes the different strategies and technologies for gene targeting.



We generated factor *FVIII* knockout mice by TALEN technology with the efficiency of 2.56-3.57% (Table 2).


Table 2. Generating factor FVIII knockout mice with TALEN genomic editing

Talens pairs	Embryos injected	Newborn	Indels detection	Mutant mouse / embryos transferred (%)	Founder Name
<b>L1+R2</b>	96	0	0	0	
<b>L1+R3</b>	171	28	1 (3.57%)	0.58	H32 ♂
<b>L2+R1</b>	120	20	0	0	
<b>L2+R2</b>	229	39	1 (2.56%)	0.44	Q83 ♂
<b>L2+R3</b>	171	26	0	0	

## Story 2. Generating factor *FVII* knockout animals with CRISPR-Cas9

Several strategies to improve specificity of the Cas9 system have been reported, such as the paired Cas9 nickase approach in which two gRNAs target adjacent sites on opposite DNA strands and each recruit a Cas9 nickase that nicks DNA instead of cutting both strands. This method can reduce off-target modifications at sites induced by single gRNA-guided Cas9.

Fu et al. reported that truncated gRNAs (tru-gRNAs) improved Cas9 nuclease specificity in U2OS.GFP and FT-HEK293 cells by shortening the gRNAs to 17/18 nt. They found that 5'-end nucleotides are not required for standard gRNA (std-gRNAs, 20 nt) activity and compensate for mismatches at unwanted positions along the gRNA target DNA interface, as shorter gRNAs are more sensitive to mismatches and therefore exhibit higher specificity. Here, we investigated the activity and specificity of tru-RGNs in inducing coagulation factor VII (FVII) gene mutations in murine cells and its efficiency in generating gene knockout (KO) mice by zygote RNA microinjection.



We generated factor *FVII* knockout mice by CRISPR-Cas9 technology with an efficiency as high as 78.9% (15/19) (Table 3).

Table 3 Generating *FVII* KO mice with gRNA and Cas9 mRNA co-injection

gRNAs	No. injected embryos	No. embryos transferred	No. recipients	No. newborns	Mutant alleles per mouse		Mutant mouse / total mouse tested (%)	Mutant mouse / embryos transferred (%)
					1(allele)	2(allele)		
F7-1	78	71	2	18	1(1)	0(0)	1/18(5.6)	1.4
tF7-1	98	98	3	38	1(1)	18(0)	19/34(55.9)	19.4
F7-2	110	106	4	24	2(2)	6(0)	8/22(36.4)	7.5
tF7-2	85	67	2	20	15(15)	0(0)	15/19(78.9)	22.4
F7-3	80	77	2	6	2(2)	0(0)	2/6(33.3)	2.6
tF7-3	130	115	3	20	8(8)	0(0)	8/20(40)	6.9

Detection *FVII* mutation frequencies by T7EI assay and subsequent sequencing (Table 4).





Table 4. Gene mutations in factor FVII knockout mice induced by CRISPR-Cas9



Founder	Target sequences (5'-3')	Indels
RGNs F7-1		
WT	AAAGGCGTGCCAACCTCACTCCTGGAGGAGCTTTGGCCCGGCTCTCTGGAGAGAGA	
O68	AAAGGC-----TTTGGCCCGGCTCTCTGGAGAGAGA	-24 bp
RGNs tF7-1		
WT	CTGTTTCTCAGTTTTTCATAACCCAGGAGGAAGCACATGGTGTCTACACAGGCAAAGGCGTGCCAACCTCACTCC	
I18	TTGTTTTCTCAGTTTTTCATAACCCAGGAGGAAGCACATGGTGTCTACACAGGCAAAGGCGTGCCAACCTCACTCC	C→T, C→T
RGNs F7-2		
WT	CCTACACAGGCAAAGGCGTGCCAACCTCACTCCTGGAGGAGCTTTGGCCCGGCTCTCT	
I31	CTTACACAGGCAAAGGCGTGCCAACCTCACTCCTGGAGGAGCTTTGGCCCGGCTCTCT	C→T
I34	CTTACACAGGCAAAGGCGGGCCAACCTCACTCCTGGAGGAGCTTTGGCCCGGCTCTCT	C→T, T→G
RGNs tF7-2		
WT	AGGCGTGCCAACCTCACTCCTGGAGGAGCTTTGGCCCGGCTCTCTGGAGAGAGAGTGCAAT	
H66	AGGC-----CCGGCTCTCTGGAGAGAGAGTGCAAT	-30 bp
H67	AGGCGTGCCAACCTCACTC-----TCTGGAGAGAGAGTGCAAT	-22 bp
H68	AGGCGTGCCAACCTCA-----GAGGAGCTTTGGCCCGGCTCTCTGGAGAGAGAGTGCAAT	-6 bp
H69	AGGCGTGCCAACCT-----TTGGCCCGGCTCTCTGGAGACAGTGCAAT	-16 bp
H71	AGGCTTG-----GCCCT-----CC-----TCTCTGGACAGAGTGCAAT	-25 bp
H72	AGGCGTGCCAACCTC-----TCTGGAGAGAGAGTGCAAT	-27 bp
H73	AGGCGTGCCAACCTCACTCCTGAAG-AGCTTTGGCCCGGCTCTCTGGAGAGAGAGTGCAAT	-1 bp
H74	AGGCGTGCCAACCTCACTCCTGGAGGAGTTTTGGCCCGGCTCTCTGGAGAGAGAGTGCAAT	-6 bp
H75	AGGCGTGCCAACCTCACTCC-----TTTGGACCGGTTCTCTGGATATAGTGCAAT	-16 bp
H76	AGGCGTGCCAACCTCACTCACTTTAGGAGCTTTGGCCCGGCTCTCTGGAGAGAGAGTGCAAT	+8 bp
	AATATGGA	
H78	AGGCGTGCCAACCTCACT-----TGGGAGAG---GG-TCTC-----AGTGCAAT	-27 bp
H79	AGGCGTGCAAACCTC-----TTTAGCTCTGCTC-----AC--AGTGCAAT	-29 bp
H82	AGGCGTGCCAACCTCACCT-C-----AACTTTGGACCGGCTCTCTGGATATAGTGCAAT	-13 bp
H83	AGGCGT-----	-167 bp
H84	AGGCGTGCCAACCTCACT-----CCGGTCTCAGGATAGAGAGAGTGCAAT	-27 bp
RGNs F7-3		
WT	CCTGGAGGAGCTTTGGCCCGGCTCTCTGGAGAGAGAGTGCAATGAGGAACACCTGG	
I66	CCTGGAGGA-----AGTGCAATGAGGAACACCTGG	-24 bp
I68	CCTGGAG-----AGAGAGTGCAATGAGGAACACCTGG	-24 bp
RGNs tF7-3		
WT	CCTGGAGGAGCTTTGGCCCGGCTCTCTGGAGAGAGAGAGTGCAATGAGGAACAGTGC	
H54	CCTGGAGGAGCTTTGGCCCGGCTCTCTGGAGAGAGAGAGTGCAATGACGAAGAGAGC	G→C, C→G, T→A
H59	CCTGGAGGAGCTTTGGCCCGGCTTCTCGGAAAAAAAAAGGGAAAGGAGAAACATTGC	+1 bp
H60	CCTGGAGGAGCTTTGGCCTGGCTCTCTGGAGAGAGAGTGCAATGAGGAACAGTGC	C→T
H62	CCTGGAG-----AGAGTGCAATGAGGAACAGTGC	-24 bp

## **We are on Nature, Nature Medicine and Nature Biotechnology**

---



To our knowledge, RLI is the only company providing transgenic rabbit and guinea pig platform worldwide. We are promoting our service in the world prestigious journals of **Nature**, **Nature Medicine** and **Nature Biotechnology**. Please see the banners below.



nature.com : Publications A-Z index : Browse by subject      Welcome back Lillian Wong       Cart : Lillian Wong : Logout

 For more information, please see  
[http://www.renovalife.com/transgenic\\_rabbits.php](http://www.renovalife.com/transgenic_rabbits.php) 

**nature** International weekly journal of science      Search    
▶ [Advanced search](#)



Home | News & Comment | Research | Careers & Jobs | Current Issue | Archive | Audio & Video | For Authors



nature.com : Publications A-Z index : Browse by subject       Lillian Wong : Logout :  Cart

 **Puzzled by how to translate the mouse results to human?  
Need animal bioreactors for producing expensive proteins?  
It's time to consider testing your hypothesis in  
rabbit and guinea pig models.** 

**nature  
medicine**      Search    
▶ [Advanced search](#)

Home | Current issue | News & comment | Research | Archive ▼ | Authors & referees ▼ | About the journal ▼

nature.com : Publications A-Z index : Browse by subject       Lillian Wong : Logout :  Cart

 **Puzzled by how to translate the mouse results to human?  
Need animal bioreactors for producing expensive proteins?  
It's time to consider testing your hypothesis in  
rabbit and guinea pig models.** 

**nature  
biotechnology**      Search    
▶ [Advanced search](#)

Home | Current issue | News & comment | Research | Archive ▼ | Authors & referees ▼ | About the journal ▼

### **Contact:**

Renova Life Inc. 7635 Timberlin Park Blvd, Suite 1216, Jacksonville, FL 32256, USA.

T: 904-379-6507. F: 301-576-5078. E: [info@renovalife.com](mailto:info@renovalife.com)

Web: [http://www.renovalife.com/transgenic\\_rabbits.php](http://www.renovalife.com/transgenic_rabbits.php)